

Adjusted by:

E-mail:

Adjusted for:

Telephone:

Date:

Save the guide to the desktop and open it from there.

Bambino 3 has an adjustable rear axle. To find a balance that suits you the rear axle can be adjusted into different positions. Note desired measure at "Settings". Bambino 3 has a seat height of 41 cm in the front and 39 at the rear, without cushion. To adjust the seat height you can chose other rearwheel sizes or a raising Kit with suitable castor fork from the list of options.
Max userweight: 60 kg for seat width 24-33

Save the guide to the "desktop" and open it with [Adobe Reader](#), for best function
Design your wheelchair by clicking in the circles in this guide.
Blue, underlined text is clickable to [show info](#) or image.
BB = seatwidth, HH= backrest height

Original design

Alternative options

Other options

[info](#)

1) Chassis-colour

- Blue** - Racing Blue
- Black** - Anodic Black
- Clear



Optional colour - A longer delivery time for optional colour

Optional colour

2) Seat width (BB)

- [art.no.](#) 24 cm 27 cm 30 cm 33 cm
- Clear



Art.no includes chassis, seat, backrest lock, anti-tip device

3) Backrest adjustable height

- [20-28 cm blue](#) 27-35 cm blue
- 523BB17 523BB18
- 20-28 cm black 27-35 cm black
- 523BB17BL 523BB18BL
- Clear



4) Rear axle

- [Carbon](#) 6 dgr
- 33030BB
- Rensa



[Carbon](#) 6 dgr +10 mm

33031BB

With 20" rearwheels +120 mm castors

[Other](#) rear axle or seat height

[Increase the seat height](#)

[List of options](#) art.

5) Castor and fork

- [120 mm](#)
- [art.no.](#) With 20" rear wheels you need to raise the rear axle 10 mm
- Clear



[90 mm](#)

Standard with 20" rearwheels

[art.nr.](#)

[Other castor fork](#)

[Excel guide](#) art.

6) Footrest

[Footrest B3](#) - foldable, angle adjustable

42341BB

Clear



7) Brake

- [Brake B3](#)
- 4630200
- Clear



[Brake blade](#) B3 PU tyre SET

4630234

8) Rear wheel - without pushrim

- [22" Standard](#)
- [art.no.](#)
- Clear



[20" Standard](#)

[art.no.](#)

[24" Standard](#)

[art.no.](#)

[Rear wheel size](#)

[Other rear wheels](#)

[List of options](#) art.

Tire - choose tire

- [Schwalbe Right Run](#) 20",22",24"
- Clear



[Rugged](#) 24"

Schwalbe [Marathon Plus](#) 22",24"

[PU](#) 20",22",24" (massive)

9) Pushrim

[Titanium](#)

190DD11 (x2 / wheelchair)

Clear DD = Wheel diameter

[Paragrip](#)

19051DDBT x2

[Tetragrip](#)

19052DDBTx2

[Maxgrip](#) 24" (contains latex)

1903060 x2

Accessories

Side guard



- [Side guard B3](#), 4235000
- Clear

Arm support



- [Arm support B3](#) with side guard 4401400
- Adjustable in height in 3 pos.

Seat width adaptor



- [Seat width adaptor](#)-pair 4235020
 - Seat width adaptor-left 4235021
 - Seat width adaptor-right 4235022
 - Clear
- Limits the back rest angle forward to 16 mm. Narrows the seat width 7,5mm on each side.

Lateral support



- [Lateral support B3](#) adjustable cpl 4236000 /pcs
- Left Right
- Position
- Standard Standard
- Narrow -9mm Narrow -9mm

Push bar



- [Push bar B3](#) 4233000
- Height, angle adjustable with Quick release

Head support



- [Head support B3 Blue](#) 4710100
 - Head support B3 Black 4710100BL
 - Clear
 - Position Low Position High
- Height, angle and depth adjustable

Cushion Wedge cushion



- [2,5 cm](#) 572BB34
- [5 cm](#) 575BB34
- Clear
- [Wedge cushion](#) 57710BB not width 24 cm

Attachment pelvic belt

- Attachment plate pelvic belt cpl 9900247 backrest 2 point attachment
- Attachment clamp pelvic belt pair 9900248 2-point attachment
- [Attachment pelvic belt 4-p](#) plate + clamp 9900250 4-point attachment
- Attachment pelvic belt 4-p with clamps 9900248 x 2 4-point attachment

Spoke protectors



- Clear
- [24" Panthera logo](#) 4500000
- [22" Panthera logo](#) 4500100
- [20" Panthera logo](#) 4500021
- [24" Without logo](#) 4500011
- [22" Without logo](#) 4500111
- [24" Splash](#) 4500324
- [22" Splash](#) 4500322
- [20" Splash](#) 4500320
- [24" Car rim](#) 4500124
- [22" Car rim](#) 4500122
- [20" Car rim](#) 4500120
- [24" Animal Planet](#) 4500424
- [22" Animal Planet](#) 4500422
- [20" Animal Planet](#) 4500420
- [24" Smiley](#) 4500224
- [22" Smiley](#) 4500222
- [20" Smiley](#) 4500220

* Tube protection front; Can be too long if footrest is in its high position. Fits when rear axle is at std balance and footrest at low pos. 16 cm or lower, high pos. 28 cm or lower, from floor to top of footplate.



Other

- [Tube protection front](#) (chassis tube) * 9900251 (2 pcs)
- [Tube protection rear](#) (backrest tube) 9900251 (1 pcs) Not with pushbar
- [Backrest wedge](#) 4710012 _____ pcs
- [Tray-table](#) 6810000
- [Magnetic holder](#) on rearwheel for tray table 6810010
- [Torso strap](#) 4380238
- [Freewheel](#) 2600102
- Freewheel-kit assembly parts 42342BB

Settings

Backrest angle



Note the measure (A) for desired backrest angle. 14-28mm. Standard=18 mm

A= _____ mm

Rear axle adjustment



Note the measure (A) for desired balance. Standard=77 mm

A= _____ mm

Footrest position



Note the height measure from the floor, (A) or (B) for the footrest.

[Bambino footplate with Feewheel - info](#)

A= _____ mm B= _____ mm

At 22" rear wheels

Measure high pos. 210-310 mm
Measure low pos. 150-220 mm

Note depth-position:

- Std position
- Moved forward position 30 mm
- High position
- Low position

Transport marking

- YES do mark the wheelchair with attachment points for transport in disabled service Art.no.9000152

NOTE! See the manual for more information.

Extra