

Adjusted by:

E-mail:

Adjusted for:

Tel:

Date:

Total weight (g):

The weight changes based on the choices made below.

U3 has a adjustable rear axle in carbon fibre. To find a balance that suits you the rear axle can be adjusted into different positions. Note desired measure at "Settings". S3 has a seat height of 47 cm in the front and 43 at the rear, without cushion. To adjust the seat height you can chose other rearwheel sizes or a raising Kit with suitable castor fork from the List of options.
Max userweight: 100 kg for seat width 33-42 and 150kg for seat width 45 cm

Save the guide to the "desktop" and open it with [Adobe Reader](#), for best function Design your wheelchair by clicking in the circles in this guide. Blue, underlined text is clickable to [show info](#) or image. **BB** = seatwidth, **HH**= backrest height

Original design

Alternative options

Other options

[info](#)

1) Chassis-colour

Black - Anodic Black

Optional colour - A higher price and

longer delivery time for optional colour, Fili in art.no. in the square



2) Seat width (BB)

[art.no.](#)

Clear

33 cm 36 cm 39 cm 42 cm 45 cm

Art.no include: chassis, seat, calf strap, backrest lock



3) Back rest

Standard
523BBHH VER2

Clear

Height (HH)

20 cm 35 cm
 25 cm 40 cm
 30 cm 45 cm

Clear

Tapered
525BBHH VER2

Curved back
524BBHH VER2

Other back rest

[List of options](#) art. _____



4) Rear axle

Carbon standard
33023BB

Clear

Carbon +15 mm
33023BC*

* seat width BB +1

Use this rear axle when chosing arm support at Accessories and/or increased seat height

Other rear axle or seat height

[Increase the seat height](#)

[List of options](#) art. _____




5) Castor and fork

90 mm
[art.no.](#)

Clear

120 mm
[art.no.](#)



6) Footrest

Footrest Titanium
36915BB

Clear

Extended, reinforced
36006BB (w. 33-36) 36006BBB (w. 39-45)

Plastic plate for footrest
36010BB

Other footrest

[List of options](#) art. _____



7) Brake

Brake S3
4630000

Clear

One hand brake
[Art.no.](#)

Other brake

[List of options](#) art. _____



8) Rear wheel - without pushrim

24" standard
[art.no.](#)

Clear

24" Spinergy X black spokes [Rear wheel sizes](#)
[art.no.](#)

25" Spinergy X black spokes
[art.no.](#)

26" Spinergy X black spokes
[art.no.](#)

Other rearwheels

[List of options](#) art. _____



Tire - choose tire

Schwalbe Right Run 24",25",26"

Clear

Rugged 24"

Schwalbe **Marathon Plus** 24",25",26"

PU 24" (massive)



9) Push rim

Titanium
190DD11 (x2 / wheelchair)

Clear

DD = Wheeldiameter

Aluminum
1904404 x2 (04=24", 05=25", 06=26")

Tetragrip
19052DDBTx2

Paragrip
19051DDBTx2

Maxgrip 22", 24"(contains latex)
1903060 x2 (60=24", 80=22")

Accessories



Side guards

- [Side guards S3](#) [Video](#)
art.no.
- Clear



Arm support

- [Arm support S3](#) [Video](#)
art.no.



Side guard Carbon

- [Side guard Carbon](#)
art.no.



Push handles

- [Height adjustable S3](#) [Video](#)
4252060 VER2
- [Height adjustable S3 +10](#)
4252063 VER2
Max. userweight 100kg
- Clear



Push handle Foldable

- [Foldable S3](#) [Video](#)
4245000



Push handle Pushbar

- [Pushbar S3](#)
42400BB VER2
Max. userweight 100kg



Anti-tip

- [Anti-tip](#) std seat height
4353000 [Video](#)

[Other seat height](#)




Cushion

- [2,5 cm](#)
572BB40
- [5 cm](#)
575BB40
- Clear

Wedge cushion

- [Wedge cushion](#)
57710BB



Spoke protectors

- [24" logo](#) [22" logo](#)
4500000 4500100
- [24" no logo](#) [22" no logo](#)
4500011 4500111
- Clear



Other

- [Tube protection front](#)
9900251 (2 pcs)
- [Tube protection rear](#)
9900251 (1 pcs)
- [Crutch holder S3](#) _____ pcs
- [9900301](#) (x1/crutch)
- [Backrest wedge](#) _____ pcs
- 4710012

Chose art.no. for Freewheel according to the **footrest height** of the wheelchair
Read more about Freewheel on the [homepage](#)

- [Freewheel](#)

- [Table S3](#)
6810020



Set-up adjustments



Back rest angle

Note measure (A) for desired back rest angle.
14-28 mm. Standard=20 mm

A= _____ mm



Rear axle adjustment

Note measure (A) for desired balance.
Standard=60 mm, at raised rear axle= 51 mm

A= _____ mm [Info - 26" rearwheel together with the S3 brake](#)



Footrest height

Note measure (A) for the footrest.

A= _____ mm



Transportation marking

- YES do mark the wheelchair with attachment points for transport in disabled service vehicle
Art.no.9000152

NOTE! See the [manual](#) for more info.

Extra